

BIAC New User Orientation

Jenn Graves | Duke-UNC Brain Imaging and Analysis Center

Agenda

- Introductions
- Facilities
- Accounts and Access
- Scheduling and Calendar Entries
- Protocol Planning
- Participant Screening and MR Safety
- Day of Scan Expectations
- BIAC Resources and Billing System
- FAQs
- Next Steps

Introductions

Welcome to BIAC!

- The Brain Imaging and Analysis Center (BIAC) was created in 1998 as the key component of the campus-wide neuroimaging initiative.
- We are now home to more than sixty researchers and staff members from throughout Duke University and UNC-Chapel Hill, and across diverse disciplines



Introductions: Dr. Todd Harshbarger & MR Technologists




Dr. Todd Harshbarger
Director of MR Operations



Jenn Graves
*Chief MR
Technologist*



Lamont Conyers




Susan Music


Facilities

Where is BIAC located?


- The main BIAC administrative offices are housed in the Davison Building, fourth floor (Duke Hospital South Green Zone)
- Duke Hospital North houses our research-dedicated MRI suites and our testing rooms.



Location within Duke Hospital




Route
to the
Scanners
from
each
lobby






MRI Entrance

Overview of the Scanner Hallway



Key Places

- 1) MRI Department Entrance
- 2) Restroom
- 3) Restroom
- 4) Lockers
- 5) MOCK I
- 6) DHNI



4) Lockers




6)MOCKI




6) DHNI



Overview of the Scanner Hallway





Key Places

7) BIAC 6

8) TEST 1

9) Restroom

10) BIAC 5

7) BIAC 6



BIAC 6 Suite



BIAC 6 Scanner



8) TEST I



10) BIAC 5

(Hallway
Entrance)



10) BIAC 5 Suite



BIAC 5 Scanner



DAV I *



*not located in the MRI hallway

Facilities Overview

- Lockers
- Test Rooms
 - MOCK I
 - DHNI
 - TEST I
 - DAV I*
- Scanners
 - BIAC 5
 - BIAC 6

Accounts and Access

To request access...

- 1) ****Must have a DHE account for BIAC access to be granted****
- 2) Email biac-newuser@duke.edu with your general access request and include:
 - a) Email
 - b) NetID
- 3) BIAC will then create your record in REDCap
 - a) *If you need Level I Access, a New User Orientation will also be set up at this time
- 4) Once your record is created an email will be sent from biac-newuser@duke.edu to the one you previously provided
 - a) Click on “BIAC User Profile” or the link to begin the survey

4) Login in and complete the survey using your NetID and email

5) If you requested Level I Access, an additional email will be sent titled, “Level I Access Survey”

a) ****DO NOT begin this survey until you have:**


- i. Scheduled a New User Orientation
- ii. Completed the MRI Safety Quiz
- iii. Completed the MRI Level I Access Signature Sheet
- iv. Completed the MRI Subject Screening form

b) This survey is where you will upload both documents and confirm you have completed the safety quiz

- a) The safety tutorials, safety quiz and both documents can be found here: [Safety | The Brain Imaging and Analysis Center \(duke.edu\)](https://www.biac.duke.edu/research/safety) or <https://www.biac.duke.edu/research/safety>

c) *the login for this survey is your NetID and email

- a) You may return to this survey, you will just need to log back in each time

 **Survey Login** ✕

Survey title: "BIAC User Profile"

Before beginning or continuing this survey, you must first log in by successfully entering the correct values below. **You must successfully enter a value for ALL the fields below.** Please note that the login is *not* case sensitive.

Net ID	<input type="text" value="test1"/>
	<input checked="" type="checkbox"/> Show value
User email	<input type="text" value="test1@duke.edu"/>
	<input checked="" type="checkbox"/> Show value

- 6) Once you have completed the “BIAC User Profile” survey, the “Level I Required Training and Documents” survey (only if need Level I access), and completed the New User Orientation – your sponsor will be contacted to confirm your access request
- 7) Once your sponsor confirms, access will then be provided

**Once you have gained access, please email biac-help@duke.edu with any questions or to fix any problems with access

Access FAQs

- Q: What is Level I Access? Do I need it?
- A: This access is for those who plan to bring participants to the scanner suites and/or use the test rooms. For Level I Access you will need to complete the New User Orientation, the MRI Safety quiz, the MRI Subject Screening Form, and the MRI Level I Access Signature Sheet.
- Q: How do I schedule a New User Orientation?
- A: Email biac-newuser@duke.edu and request a New User Orientation session. All sessions are currently taught over Zoom for approximately, 1.5 hours.
- Q: I already have data access but need to get scheduling access, who do I email?
- A: Email biac-help@duke.edu with your access request with your sponsor Cced so your request can be confirmed.

FAQs Continued

- Q: Can I access the REDCap surveys from anywhere or do I have to be on the Duke network?
- A: The BIAC User Profile survey and the Level I Access survey through REDCap can be taken anywhere. You will just need to login into the quiz using the NetID and email you previously provided
- Q: I have now have access to the BIAC Scanner Schedule, but I cannot login to view it...
- A: Make sure you are on the DHE network or on a DHE computer. If you are working remotely, you may need to request DHE access.

DHE Access

- Not the same as your Duke Account
- To get DHE Access
 - Need to have someone with a DHE appointment (faculty/staff) sponsor you at this site: [Sponsored Accounts \(duke.edu\)](#)
- ****MUST HAVE a DHE account to gain BIAC access****

The screenshot shows the 'Duke Sponsored Accounts' website. At the top, there is a dark blue header with the text 'Duke Sponsored Accounts' in white and yellow. Below the header, there are two navigation tabs: 'Request/Manage' (which is active) and 'Policies'. The main content area is titled 'Request and Manage Accounts' and is divided into three columns:

- DHE Access:** For Duke staff, faculty, or students who need to be sponsored for Duke Health Enterprise access. A dark blue button reads: 'Request a Duke Health DHE Account for a Duke employee or Duke student *'.
- Guest (Affiliate) Accounts:** For non-Duke users needing physical or electronic access at Duke. A dark blue button reads: 'Request a new guest account for a non-Duke user'.
- Non-person Accounts:** For test and service accounts to support Duke services. A white button reads: 'Manage my test and service accounts (0)'. Below this, a message states: 'You do not have access to request a non-person account. More information [here](#).'

DHE Account Request Form

Request a Duke Health DHE Account for a Duke employee or Duke student

1. Select DHE account type*



- Collaborator
- Contractor
- Duke Health Employee
- Student, non-SoM/SoN
- Student, SoM/SoN
- Training
- Volunteer
- Other

2. Set Account Preferences



User Information

NetID of user needing Duke Health DHE Account*:

Business justification*:

User is a minor (less than 18 years old)

Account Length

Account(s) will expire on **04/29/2024** unless renewed.

[Or, choose an earlier expiration date.](#)

I agree to the [terms of DHE account sponsorship.](#)

Submit

Scheduling a Scan

Step-by-step

Scheduling

- Most often, the experimenter (i.e. the person running the participant) schedules a session on a BIAC scanner
- It is the expectation that you have a participant when you reserve a slot
- Last minute scheduling
 - If scheduling a scan less than 24 hours before their start time or on the same day, please email BIAC-Schedule@duke.edu to ensure coverage will be provided

How far in advance can you schedule?

Experiment Information for Demo.01

Go back to [main calendar view](#).
Show [all experiments](#) user "jel76" has access to.

Summary

This experiment is currently valid.
This experiment's scheduling window allows you to schedule through 7/5/2023.
This experiment has 269 incomplete calendar entries ([see below](#)).

You may contact the [BIAC Scientific Review Committee](#) to request more time or to [renew your experiment](#).
See [SOP109: Approval for BIAC Scanner Use](#) for more details.

Experiment	Scanner	Investigator	Responsible Investigator	Department	Start Date	End Date	
Demo.01		Allen Song	Allen Song	BIAC	12/1/2000	1/1/2024	Demonstrations using MR facilities

Experiment	Time Allocated (hours)	Time Left (hours)	Time Used (hours)	Lost Time [Claimed] (hours)	Total Time (hours)	Scheduling Window (5/3/2023 - 5/24/2023)			Scheduling Window (5/25/2023 - 7/5/2023)		
						Window Size (days)	Max Study Size (hours)	Max Weekly Time (hours)	Window Size (days)	Max Study Size (hours)	Max Weekly Time (hours)
Demo.01	1000	190.5	809.5	8.5 [8.5]	818	21	4	None	42	4	6

- 2 windows:

- 21 days out


- No limit on number of weekly hours, but can only schedule 4 hours at a time

- 42 days out

- Weekly limit of 4 hours (can be increased if needed)
 - Can schedule 4 hours at a time

Scheduling Step 1: Locate the calendar page

- Go to <https://www.biac.duke.edu/>
- Click on the “Services” tab



Scheduling Step 2: Access the Calendar

- Once in the Service tab, click on the Calendar icon
- *Must be on the Duke Hospital (DHE) network or computer to access the calendar

The screenshot shows a web browser window with the URL biac.duke.edu/services. The page header includes the Duke-UNC Brain Imaging and Analysis Center logo and name. A navigation menu contains links for About Us, Research, Education, People, Facilities, and Services. Below the menu, there are three service tiles: Scanner Calendar, Trouble Ticket, and Exam Tracker. A red arrow points to the Scanner Calendar tile.

BIAC Scanner Schedule

BIAC Scanner Schedule

[Experiment Info Page](#)
[MR Tech Schedule](#)

Scanner:

Start Time:

Stop Time:

Update View

March

S	M	T	W	T	F	S
			1	2	3	4
5	6	7	8	9	10	11
12	13	14	15	16	17	18
19	20	21	22	23	24	25
26	27	28	29	30	31	

Experimenters must complete all calendar entries for their experiments.

Color Key	
Info Complete	Canceled
Info Missing	Overlap

Calendar version 1.7.15 [2022.03.16]
 by [Francis Favorini](#).

	Sunday 3/19/2023	Monday 3/20/2023	Tuesday 3/21/2023	Wednesday 3/22/2023	Thursday 3/23/2023	Friday 3/24/2023	Saturday 3/25/2023
	Add item	Add item	Add item	Add item	Add item	normal tech coverage 8-3:30 single tech coverage 3:30-5	no tech
8:00 AM							
9:00 AM						Null.01	
10:00 AM		Null.01 Demo.01 Gloria Meiers			CHI.01 Ryan Bell	NOGGIN.01 Jeffrey Browndyke 9:00 AM-10:30 AM	
11:00 AM			ALADDIN.01 Jeffrey Browndyke	Null.01 UserTest.01 Gloria Meiers	Null.01 Ryan Bell	Null.01	NCANDA.04 Casey Pearce 11:00 AM-12:30 PM
12:00 PM		NoTech.01 None 12:00 PM-2:00 PM 00000		QA.04 Jeffrey Browndyke	NOGGIN.01 Jeffrey Browndyke 11:30 AM-1:00 PM		
1:00 PM			Null.01	QA.04	Null.01	NoTech.01	
2:00 PM			NoTech.01	NoTech.01	Null.01	None	
3:00 PM			QA.04	QA.04	Null.01	Null.01	
4:00 PM			QA.04	Cleaning.01	CHI.01 Ryan Bell	ADRC.01 Emily Maingi	
5:00 PM			Null.01 Ryan Bell	UserTest.01 Trong-Kha Truong		Null.01	
6:00 PM			CHI.01	SysDevelop.53		Emily Maingi	
7:00 PM			QA.04				
8:00 PM		Cleaning.01 SysDevelop.53					
9:00 PM			Cleaning.01	SysDevelop.53			

BIAC Scanner Schedule

2 BIAC Scanner Schedule

[Experiment Info Page](#)
[MR Tech Schedule](#)

Scanner:
 Start Time:
 Stop Time:

Update View

S	M	T	W	T	F	S
		1	2	3	4	
5	6	7	8	9	10	11
12	13	14	15	16	17	18
19	20	21	22	23	24	25
26	27	28	29	30	31	

Experimenters must complete all calendar entries for their experiments.

Color Key	
Info Complete	Canceled
Info Missing	Overlap

Calendar version 1.7.15 [2022.03.16] by Francis Favorini.

	Sunday 3/19/2023	Monday 3/20/2023	Tuesday 3/21/2023	Wednesday 3/22/2023	Thursday 3/23/2023	Friday 3/24/2023	Saturday 3/25/2023
	Add item	Add item	Add item	Add item	Add item	normal tech coverage 8-3:30 single tech coverage 3:30-5	no tech
8:00 AM	1						
9:00 AM						Null.01	
10:00 AM		Null.01 Demo.01 Gloria Meiers			Null.01 CHI.01 Ryan Bell	NOGGIN.01 Jeffrey Browndyke 9:00 AM-10:30 AM	
11:00 AM			ALADDIN.01 Jeffrey Browndyke	Null.01 UserTest.01 Gloria Meiers	Null.01 NOGGIN.01 Jeffrey Browndyke 11:30 AM-1:00 PM	Null.01 NCANDA.04 Casey Pearce 11:00 AM-12:30 PM	
12:00 PM		NoTech.01 None 12:00 PM-2:00 PM 00000	Null.01	QA.04 QA.04	Null.01 NoTech.01 Jeffrey Browndyke 11:30 AM-1:00 PM		
1:00 PM			Null.01	Null.01	Null.01	NoTech.01 None	
2:00 PM			NoTech.01	NoTech.01 Null.01	Null.01		
3:00 PM			QA.04 QA.04	Null.01	Null.01	Null.01	
4:00 PM			Null.01 CHI.01 Ryan Bell	Cleaning.01 UserTest.01 Trong-Kha Truong	CHI.01 Ryan Bell	ADRC.01 Emily Maingi	
5:00 PM			CHI.01 Ryan Bell	SysDevelop.53		Null.01 Emily Maingi	
6:00 PM			QA.04				
7:00 PM							
8:00 PM		Cleaning.01 SysDevelop.53					
9:00 PM			Cleaning.01 SysDevelop.53				

Scheduling Step 3: Select the scanner and choose an open time slot

BIAC Scanner Schedule

1

[Experiment Info Page](#)
[MR Tech Schedule](#)

Scanner: **BIAC5** ▼

Start Time: 8:00 AM
 Stop Time: 10:00 PM

Update View

3/23/2023 Today

March 2023


S	M	T	W	T	F	S
		1	2	3	4	
5	6	7	8	9	10	11
12	13	14	15	16	17	18
19	20	21	22	23	24	25
26	27	28	29	30	31	

Experimenters must complete all calendar entries for their experiments.

Color Key	
Info Complete	Canceled
Info Missing	Overlap

Calendar version 1.7.15 [2022.03.16]
 by Francis Favorini.

	Sunday 3/19/2023	Monday 3/20/2023	Tuesday 3/21/2023	Wednesday 3/22/2023	Thursday 3/23/2023	Friday 3/24/2023	Saturday 3/25/2023
	Add item	Add item	Add item	Add item	Add item	normal tech coverage 8-3:30 single tech coverage 3:30-5	no tech
8:00 AM						Null.01	
9:00 AM						NOGGIN.01	
10:00 AM		Null.01 Demo.01 Gloria Meiers	Null.01	Null.01	CHI.01 Ryan Bell	Jeffrey Browndyke 9:00 AM-10:30 AM Null.01	
11:00 AM			ALADDIN.01 Jeffrey Browndyke	UserTest.01 Gloria Meiers	Null.01 NOGGIN.01 Jeffrey Browndyke 11:30 AM-1:00 PM Canceled	NCANDA.04 Casey Pearce 11:00 AM-12:30 PM	
12:00 PM		NoTech.01 None 12:00 PM-2:00 PM 00000	Null.01	QA.04 QA.04	Null.01	NoTech.01	
1:00 PM			Null.01	Null.01	Null.01	None	
2:00 PM			NoTech.01	Null.01	Null.01		
3:00 PM			QA.04 QA.04	Cleaning.01 UserTest.01 Trong-Kha Truong	Null.01 CHI.01 Ryan Bell	Null.01	
4:00 PM			Null.01 CHI.01 Ryan Bell	SysDevelop.53		ADRC.01 Emily Maingi Null.01 Emily Maingi	
5:00 PM			QA.04				
6:00 PM							
7:00 PM							
8:00 PM		Cleaning.01 SysDevelop.53					
9:00 PM			Cleaning.01 SysDevelop.53				

2

 Free Slot

Scheduling Step 4: Enter Study details

Study Information

biac-web.dhe.duke.edu/calendar/recdetail.asp?Date=3/23/2023&ScheduledStartTime=6%3A30%20PM&Scanner=BIAC5

Study Details

[Experiment Info Page](#)

Date: 3/23/2023 Thursday

*Scanner: BIAC5

Subject: [Redacted]

Notes: [Redacted]

*Start Time: 6:30 PM

*Experiment: [Redacted]

Change Subject

*End Time: [Redacted]

*Experimenter: [Redacted]

*Subject Type: None

*Required fields.

Add Don't Add

1

2

3

4

1) Fill out start and end time of scan

2) Experiment

3) Experimenter

4) Subject type

Scheduling Step 5: Enter Study details – Patient Information

The screenshot displays a web application interface for entering study details. The main window, titled "Study Information", shows a form with the following fields:

- Date: 3/23/2023 (Thursday)
- *Scanner: BIACS
- *Start Time: 6:30 PM
- *Experiment: (empty dropdown)
- Subject: (empty text field)
- Notes: (empty text area)

Below the form, there are buttons for "Add" and "Don't Add". A "Change Subject" button is located below the "Subject" field. A red arrow points from this button to a secondary window titled "Subject Information - Google Chrome".

The "Subject Information" window contains the following fields:

- *Last Name: doe
- *Birth Date: 1/1/2001
- *Required fields. (dropdown menu)
- Buttons: Find, No Subject

Scheduling Step 6: Save

Study Information x +

biac-web.dhe.duke.edu/calendar/recdetail.asp?Date=3/23/2023&ScheduledStartTime=6%3A30%20PM&Scanner=BIAC5

Study Details

[Experiment Info Page for Demo.01](#)

Date: 3/23/2023 Thursday

*Scanner: BIAC5

Subject: [Redacted]

*Start Time: 6:30 PM

*Experiment: Demo.01

*End Time: 7:00 PM

*Experimenter: None

*Subject Type: Human

Change Subject

Notes: [Empty text area]

required fields.

Add Don't Add

Once all the required information is selected from the drop downs, click “Add” to save the scan on the schedule

Scheduling Step 7: Make sure the scan is on the schedule

BIAC Scanner Schedule

biac-web.dhe.duke.edu/calendar/default.asp

BIAC Scanner Schedule

[Experiment Info Page](#)
[MR Tech Schedule](#)

Scanner:
 Start Time:
 Stop Time:

S	M	T	W	T	F	S
		1	2	3	4	
5	6	7	8	9	10	11
12	13	14	15	16	17	18
19	20	21	22	23	24	25
26	27	28	29	30	31	

Experimenters must complete all calendar entries for their experiments.

Color Key
 Info Complete Canceled
 Info Missing Overlap

Calendar version 1.7.15 [2022.03.16] by [Francis Favorini](#).

	Sunday 3/19/2023	Monday 3/20/2023	Tuesday 3/21/2023	Wednesday 3/22/2023	Thursday 3/23/2023	Friday 3/24/2023	Saturday 3/25/2023
	Add item	Add item	Add item	Add item	Add item	normal tech coverage 8-3:30 single tech coverage 3:30-5	no tech
8:00 AM						Null.01	
9:00 AM						NOGGIN.01 Jeffrey Browndyke 9:00 AM-10:30 AM	
10:00 AM		Null.01 Demo.01 Gloria Meiers	Null.01	Null.01	CHI.01 Ryan Bell	Null.01	
11:00 AM			ALADDIN.01 Jeffrey Browndyke	UserTest.01 Gloria Meiers	Null.01	NCANDA.04 Casey Pearce 11:00 AM-12:30 PM	
12:00 PM		NoTech.01 None 12:00 PM-2:00 PM 00000	Null.01	QA.04 QA.04	NOGGIN.01 Jeffrey Browndyke 11:30 AM-1:00 PM Canceled		
1:00 PM			Null.01	Null.01	Null.01	NoTech.01 None	
2:00 PM			NoTech.01	NoTech.01	Null.01		
3:00 PM			QA.04		Null.01	Null.01	
4:00 PM			QA.04	Cleaning.01	CHI.01 Ryan Bell	ADRC.01 Emily Maingi	
5:00 PM			Null.01 CHI.01 Ryan Bell	UserTest.01 Trong-Kha Truong SysDevelop.53		Null.01 Emily Maingi	
6:00 PM			QA.04				
7:00 PM					Demo.01		
8:00 PM		Cleaning.01 SysDevelop.53					
9:00 PM			Cleaning.01				
			SysDevelop.53				

What to do if you can't schedule your scan

Locate the Experiment Info page

BIAC Scanner Schedule

[Experiment Info Page](#)
[MR Tech Schedule](#)

Scanner:

Start Time:

Stop Time:

Study Information x +

← → ↻ 🏠 biac-web.dhe.duke.edu/calendar/recdetail.asp?Date=3/27/2023&ScheduledStartTime=4%3A30%20PM&Scanner=BIAC5

Study Details

[Experiment Info Page for Demo.01](#)

Date	<input type="text" value="3/27/2023"/>	<input type="text" value="Monday"/>	*Start Time	<input type="text" value="4:30 PM"/>
*Scanner	<input type="text" value="BIAC5"/>		*Experiment	<input type="text" value="Demo.01"/>
Subject	<input type="text"/>		<input type="button" value="Change Subject"/>	

Experiment Information Page

Experiment Information for Demo.01

Go back to [main calendar view](#).
Show [all experiments](#) user "jel76" has access to.

Summary

This experiment is currently valid.
This experiment's scheduling window allows you to schedule through 7/5/2023.
This experiment has 269 incomplete calendar entries ([see below](#)).

You may contact the [BIAC Scientific Review Committee](#) to request more time or to [renew your experiment](#).
See [SOP109: Approval for BIAC Scanner Use](#) for more details.

Experiment	Scanner	Investigator	Responsible Investigator	Department	Start Date	End Date	Scheduling Window (5/3/2023 - 5/24/2023)			Scheduling Window (5/25/2023 - 7/5/2023)		
Demo.01		Allen Song	Allen Song	BIAC	12/1/2000	1/1/2024	Window Size (days)	Max Study Size (hours)	Max Weekly Time (hours)	Window Size (days)	Max Study Size (hours)	Max Weekly Time (hours)
Demo.01							21	4	None	42	4	6

Experiment	IRB Protocol	Principal Investigator	Start Date	End Date
Demo.01	0000-00-3	Allen Song	1/1/1999	1/1/2100
Demo.01	0000-00-4	Allen Song	1/1/1999	1/1/2100

- 1. Has your experiment expired?
- 2. Has the experiment run out of hours?
- 3. Are you outside of the scheduling window?
- 4. Is the IRB protocol expired?

IRB Renewal: Calendar Entries

- Incomplete calendar entries will cause your BIAC experiment to be ineligible for renewal
- Experiment info page – will show all incomplete entries
 - This page also provides the link to the steps to take to renew your experiment when the time comes
 - Link:
<https://www.biac.duke.edu/research/scheduling-research-biac>

Please fix the following incomplete calendar entries!

Experiment	Date	Start	End	Scanner	Experimenter
Demo.01	1/27/2023	4:30 PM	5:30 PM	BIAC5	Todd Harshbarger
Demo.01	1/5/2023	2:30 PM	5:30 PM	DHN1	Rachael Wright
Demo.01	11/15/2022	3:30 PM	4:30 PM	BIAC5	Ryan Bell

Experiment Information for Demo.01

Go back to [main calendar view](#).
Show [all experiments](#) user "jel76" has access to.

Summary

This experiment is currently valid.
This experiment's scheduling window allows you to schedule through 7/6/2023.
This experiment has 269 incomplete calendar entries ([see below](#)).

You may contact the [BIAC Scientific Review Committee](#) to request more time or to [renew your experiment](#).
See [SOP109: Approval for BIAC Scanner Use](#) for more details.

Scheduling a Null

Step-by-step

Null.01

- When to use:
 - To avoid double booking when single tech coverage
 - To schedule the required 30-minute cleaning session before the scan
- When not to use:
 - To book a Usertest
 - To book a paid scan
 - When multiple techs are available

***In the event of a double-booking, the study that scheduled first (but failed to add a Null) will be rescheduled**

Null.01

- Scheduled the same as a scan
- DOES NOT contain participant information
- If you delete or cancel a scan, remove the Nulls as well

To schedule a Null for cleaning...

BIAC Scanner Schedule

Experiment Info Page
MR Tech Schedule

Scanner: BIAC5
Start Time: 8:00 AM
Stop Time: 10:00 PM

Update View

3/23/2023 Today
March 2023

S	M	T	W	T	F	S
			1	2	3	4
5	6	7	8	9	10	11
12	13	14	15	16	17	18
19	20	21	22	23	24	25
26	27	28	29	30	31	

Color Key
Info Complete Canceled
Info Missing Overlap

Calendar version 1.7.15 [2022.03.16]
by Francis Favonni.

	Sunday 3/19/2023	Monday 3/20/2023	Tuesday 3/21/2023	Wednesday 3/22/2023	Thursday 3/23/2023	Friday 3/24/2023	Saturday 3/25/2023
	Add item	Add item	Add item	Add item	Add item	normal tech coverage 8-3:30 single tech coverage 3:30-5	no tech
8:00 AM						Null.01	
9:00 AM		Null.01			Null.01	NOGGIN.01 Jeffrey Browndyke 9:00 AM-10:30 AM	
10:00 AM		Demo.01 Gloria Meiers	Null.01	Null.01	CHI.01 Ryan Bell	Null.01	
11:00 AM			ALADDIN.01 Jeffrey Browndyke	UserTest.01 Gloria Meiers	Null.01	NCANDA.04 Casey Pearce 11:00 AM-12:30 PM	
12:00 PM		NoTech.01 None 12:00 PM-2:00 PM 00000	Null.01	QA.04 QA.04	NOGGIN.01 Jeffrey Browndyke 11:30 AM-1:00 PM Canceled		
1:00 PM			Null.01	Null.01	Null.01	NoTech.01 None	
2:00 PM			Null.01	NoTech.01	NoTech.01		
3:00 PM			QA.04	QA.04	Null.01	Null.01	
4:00 PM			QA.04	Cleaning.01	CHI.01 Ryan Bell	ADRC.01 Emily Maingi	
5:00 PM			Null.01	UserTest.01 Trong-Kha Truong		Null.01	
6:00 PM			CHI.01 Ryan Bell	SysDevelop.53		Null.01 Emily Maingi	
7:00 PM			QA.04		Demo.01		
8:00 PM		Cleaning.01 SysDevelop.53					
9:00 PM			Cleaning.01				

Free Slot

Experimenters must complete all calendar entries for their experiments.

Click on the time slot right before your scan

Study Information x +

biac-web.dhe.duke.edu/calendar/recdetail.asp?Date=3/23/2023&ScheduledStartTime=6%3A00%20PM&Scanner=BIAC5

Study Details

[Experiment Info Page for Null.01](#)

Date	3/23/2023 Thursday	*Start Time	6:00 PM	*End Time	6:30 PM
*Scanner	BIAC5	*Experiment	Null.01	*Experimenter	None
Subject		<input type="button" value="Change Subject"/>		*Subject Type	None
Notes	<input type="text"/>				

*Required fields.

- 1) End time will be 30 minutes after your start time
- 2) Select “Null.01” for the Experiment
- 3) Select your name for the Experimenter
- 4) Click “Add” to save

*Nulls are not for scheduling participants so “Change Subject” and Subject Type” options will NOT be available

BIAC Scanner Schedule

[Experiment Info Page](#)
[MR Tech Schedule](#)

Scanner:
 Start Time:
 Stop Time:

S	M	T	W	T	F	S
			1	2	3	4
5	6	7	8	9	10	11
12	13	14	15	16	17	18
19	20	21	22	23	24	25
26	27	28	29	30	31	

Experimenters must complete all calendar entries for their experiments.

Color Key	
Info Complete	Canceled
Info Missing	Overlap

Calendar version 1.7.15 [2022.03.16]
 by [Francis Favorini](#).

	Sunday 3/19/2023	Monday 3/20/2023	Tuesday 3/21/2023	Wednesday 3/22/2023	Thursday 3/23/2023	Friday 3/24/2023	Saturday 3/25/2023
	Add item	Add item	Add item	Add item	Add item	normal tech coverage 8-3:30 single tech coverage 3:30-5	no tech
8:00 AM							
9:00 AM						Null.01	
10:00 AM		Null.01 Demo.01 Gloria Meiers			CHI.01 Ryan Bell	NOGGIN.01 Jeffrey Browndyke 9:00 AM-10:30 AM	
11:00 AM			Null.01	Null.01	Null.01	Null.01	
12:00 PM		NoTech.01 None 12:00 PM-2:00 PM 00000	ALADDIN.01 Jeffrey Browndyke	UserTest.01 Gloria Meiers	NOGGIN.01 Jeffrey Browndyke 11:30 AM-1:00 PM Canceled	NCANDA.04 Casey Pearce 11:00 AM-12:30 PM	
1:00 PM			Null.01	QA.04	Null.01	NoTech.01	
2:00 PM			Null.01	NoTech.01	NoTech.01	None	
3:00 PM			NoTech.01	QA.04	Null.01	Null.01	
4:00 PM			QA.04	Cleaning.01	CHI.01 Ryan Bell	ADRC.01 Emily Maingi	
5:00 PM			QA.04	UserTest.01 Trong-Kha Truong		Null.01 Emily Maingi	
6:00 PM			CHI.01 Ryan Bell	SysDevelop.53			
7:00 PM			QA.04		Null.01 Demo.01		
8:00 PM		Cleaning.01					
9:00 PM		SysDevelop.53	Cleaning.01	SysDevelop.53			

*Null.01 slots are only added to the slots BEFORE the scans, you do not need to add one after the scan

To schedule a Null during single tech coverage...

BIAC Scanner Schedule

biac-web.dhe.duke.edu/calendar/default.asp

BIAC Scanner Schedule

[Experiment Info Page](#)
[MR Tech Schedule](#)

Scanner: **BIACS**
Start Time: 8:00 AM
Stop Time: 10:00 PM

Update View

3/23/2023 Today
March 2023

S	M	T	W	T	F	S
			1	2	3	4
5	6	7	8	9	10	11
12	13	14	15	16	17	18
19	20	21	22	23	24	25
26	27	28	29	30	31	

Experimenters must complete all calendar entries for their experiments.

Color Key
Info Complete **Canceled**
Info Missing **Overlap**

Calendar version 1.7.15 [2022.03.16]
by [Francis Favorini](#).

	Sunday 3/19/2023	Monday 3/20/2023	Tuesday 3/21/2023	Wednesday 3/22/2023	Thursday 3/23/2023	Friday 3/24/2023	Saturday 3/25/2023
	Add item	Add item	Add item	Add item	Add item	normal tech coverage 8-3:30 single tech coverage 3:30-5	no tech
8:00 AM						Null.01	
9:00 AM		Null.01			Null.01	NOGGIN.01 Jeffrey Browndyke 9:00 AM-10:30 AM	
10:00 AM		Demo.01 Gloria Meiers	Null.01	Null.01	CHI.01 Ryan Bell	Null.01	
11:00 AM			ALADDIN.01 Jeffrey Browndyke	UserTest.01 Gloria Meiers	Null.01	NCANDA.04 Casey Pearce 11:00 AM-12:30 PM	
12:00 PM		NoTech.01 None 12:00 PM-2:00 PM 00000	Null.01	QA.04	NOGGIN.01 Jeffrey Browndyke 11:30 AM-1:00 PM Canceled		
1:00 PM				Null.01	Null.01	NoTech.01	
2:00 PM			NoTech.01	NoTech.01	Null.01	None	
3:00 PM			QA.04	QA.04	Null.01	Null.01	
4:00 PM			QA.04	Cleaning.01	CHI.01 Ryan Bell	ADRC.01 Emily Maingi	
5:00 PM			Null.01	UserTest.01 Trong-Kha Truong		Null.01 Emily Maingi	
6:00 PM			CHI.01 Ryan Bell	SysDevelop.53			
7:00 PM			QA.04		Null.01		
8:00 PM		Cleaning.01			Demo.01		
9:00 PM		SysDevelop.53		Cleaning.01			
			SysDevelop.53				

Check the MRI Tech Schedule to see normal coverage

MR Tech Schedule

- Normal tech coverage
- Be sure to check the notes (in blue) at the top of the scanner schedule to see exceptions to the normal tech schedule

Friday
3/24/2023

[normal tech coverage 8-3:30](#)
[single tech coverage 3:30-5](#)

BIAC Scanner Schedule

[Experiment Info Page](#)
[MR Tech Schedule](#)

Scanner:

Start Time:

Stop Time:



Susan ■
 Lamont ■
 Jennifer ■

MR Tech Schedule (Effective November 2021)

Time	Monday	Tuesday	Wednesday	Thursday	Friday	Saturday
8:00 AM - 9:00 AM						***
9:00 AM - 10:00 AM	■	■	■	■	■	■
10:00 AM - 11:00 AM	■	■	■	■	■	■
11:00 AM - 12:00 PM	■	■	■	■	■	■
12:00 PM - 1:00 PM	■	■	■	■	■	■
1:00 PM - 2:00 PM	■	■	■	■	■	■
2:00 PM - 3:00 PM	■	■	■	■	■	■
3:00 PM - 4:00 PM	■	■	■	■	■	■
4:00 PM - 5:00 PM	■	■	■	■	■	■
5:00 PM - 6:00 PM	■	■	■	■	■	■
6:00 PM - 7:00 PM	■	■	■	■	■	■
7:00 PM - 8:00 PM	■	■	■	■	■	■

*** - We have returned to alternate Saturday coverage. Check the calendar for availability. For weeks with Saturday coverage, only two techs will be available on Thursdays (shown in blue and yellow above). When there is no Saturday coverage, there will be additional Thursday coverage (shown in green).

To schedule a Null during single tech coverage...

BIAC Scanner Schedule

biac-web.dhe.duke.edu/calendar/default.asp

BIAC Scanner Schedule

[Experiment Info Page](#)
[MR Tech Schedule](#)

Scanner: **BIACS**

Start Time: 8:00 AM
Stop Time: 10:00 PM

Update View

3/23/2023 Today

March 2023

S	M	T	W	T	F	S
			1	2	3	4
5	6	7	8	9	10	11
12	13	14	15	16	17	18
19	20	21	22	23	24	25
26	27	28	29	30	31	

Experimenters must complete all calendar entries for their experiments.

Color Key

- Info Complete **Canceled**
- Info Missing **Overlap**

Calendar version 1.7.15 [2022.03.16] by Francis Favorini.

	Sunday 3/19/2023	Monday 3/20/2023	Tuesday 3/21/2023	Wednesday 3/22/2023	Thursday 3/23/2023	Friday 3/24/2023	Saturday 3/25/2023
	Add item	Add item	Add item	Add item	Add item	normal tech coverage 8-3:30 single tech coverage 3:30-5	no tech
8:00 AM						Null.01	
9:00 AM		Null.01			Null.01	NOGGIN.01 Jeffrey Browndyke 9:00 AM-10:30 AM	
10:00 AM		Demo.01 Gloria Meiers	Null.01	Null.01	CHI.01 Ryan Bell	Null.01	
11:00 AM			ALADDIN.01 Jeffrey Browndyke	UserTest.01 Gloria Meiers	Null.01	NCANDA.04 Casey Pearce 11:00 AM-12:30 PM	
12:00 PM		NoTech.01 None 12:00 PM-2:00 PM 00000	Null.01	QA.04	NOGGIN.01 Jeffrey Browndyke 11:30 AM-1:00 PM Canceled		
1:00 PM			Null.01	QA.04	Null.01	NoTech.01	
2:00 PM			NoTech.01	Null.01	Null.01	None	
3:00 PM			QA.04	QA.04	Null.01	Null.01	
4:00 PM			QA.04	Cleaning.01	CHI.01 Ryan Bell	ADRC.01 Emily Maingi	
5:00 PM			Null.01	UserTest.01 Trong-Kha Truong		Null.01 Emily Maingi	
6:00 PM			CHI.01 Ryan Bell	SysDevelop.53			
7:00 PM			QA.04		Null.01		
8:00 PM		Cleaning.01			Demo.01		
9:00 PM		SysDevelop.53					
			Cleaning.01				
			SysDevelop.53				

Book Null for both scan time and cleaning slot

To schedule a Null during single tech coverage...

BIAC Scanner Schedule

Experiment Info Page
MR Tech Schedule

Scanner:

Start Time:

Stop Time:

S	M	T	W	T	F	S
			1	2	3	4
5	6	7	8	9	10	11
12	13	14	15	16	17	18
19	20	21	22	23	24	25
26	27	28	29	30	31	

Experimenters must complete all calendar entries for their experiments.

Color Key

- Info Complete (Blue)
- Canceled (Red)
- Info Missing (Orange)
- Overlap (Purple)

Calendar version 1.7.15 [2022.03.16] by Francis Favorini.

	Sunday 3/19/2023	Monday 3/20/2023	Tuesday 3/21/2023	Wednesday 3/22/2023	Thursday 3/23/2023	Friday 3/24/2023	Saturday 3/25/2023
	Add item	Add item	Add item	Add item	Add item	normal tech coverage 8-3:30 single tech coverage 3:30-5	no tech
8:00 AM							
9:00 AM							
10:00 AM			Null.01	Null.01		Null.01	
11:00 AM			SchemRep.01 Erin Welch	SchemRep.01 Erin Welch		NOGGIN.01 Jeffrey Browndyke 10:00 AM-11:30 AM	
12:00 PM		Null.01	Null.01	Null.01	Null.01	Cleaning.01	
		IndBlast.01 Christine Baird 12:00 PM-1:30 PM	UserTest.01 John Graner	memdm.01 Simone Compton 12:00 PM-1:30 PM 04270	IntractEP.01 Prachi Parikh		
1:00 PM			Cleaning.01		Null.01	SysService.01 Christopher Petty	
2:00 PM		Null.01	NoTech.01	Null.01	Null.01	UserTest.01 Amanda Nelli	
		PROP.01 Amanda Nelli		QA.04 None	PROP.01 Amanda Nelli		
3:00 PM			QA.04 None	BrainSpine.01 Su Hyoun Park		Null.01 None 2:30 PM-5:00 PM 00000	
4:00 PM		UserTest.01 John Graner					
5:00 PM							
6:00 PM			QA.04		Null.01 None		
7:00 PM						Null.01 None 6:00 PM-7:00 PM 00000	
8:00 PM							
9:00 PM							

Book Null for both scan time and cleaning slot

Single Tech Scheduling Example...Saturday

BIAC Scanner Schedule

biac-web.dhe.duke.edu/calendar/default.asp

BIAC Scanner Schedule

[Experiment Info Page](#)
[MR Tech Schedule](#)

Scanner:
 Start Time:
 Stop Time:

S	M	T	W	T	F	S
	1	2	3	4	5	6
7	8	9	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30	31			

Experimenters must complete all calendar entries for their experiments.

Color Key
 Info Complete Canceled
 Info Missing Overlap

Calendar version 1.7.15 [2022.03.16]
by [Francis Favorini](#).

	Sunday 5/7/2023	Monday 5/8/2023	Tuesday 5/9/2023	Wednesday 5/10/2023	Thursday 5/11/2023	Friday 5/12/2023	Saturday 5/13/2023
8:00 AM	Add item	Add item	Add item	Add item	Add item	Add item	Add item
9:00 AM				Null.01			
10:00 AM				NOGGIN.01 Jeffrey Browndyke 9:30 AM-11:00 AM	BrainSpine.01 Anne Baker 10:00 AM-12:30 PM		
11:00 AM				NOGGIN.01		NoTech.01 None	BrainSpine.01 Anne Baker 11:00 AM-1:00 PM
12:00 PM		NoTech.01 None				Null.01	
1:00 PM			Null.01		SchemRep.01 Margaret McAllister	IndBlast.01 Christine Baird 12:30 PM-2:00 PM	Null.01 Victoria Swaine 1:00 PM-2:30 PM 00000
2:00 PM		BrainSpine.01 Anne Baker 2:00 PM-4:30 PM	SchemRep.01 Margaret McAllister	BrainSpine.01 Anne Baker 2:00 PM-4:30 PM	Null.01		
3:00 PM					NoTech.01		
4:00 PM					Null.01		
5:00 PM							
6:00 PM							
7:00 PM							
8:00 PM							
9:00 PM							

Single Tech Scheduling Examples...Saturday

BIAC Scanner Schedule

[Experiment Info Page](#)
[MR Tech Schedule](#)

Scanner:

Start Time:

Stop Time:

May


S	M	T	W	T	F	S
	1	2	3	4	5	6
7	8	9	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30	31			

Experimenters must complete all calendar entries for their experiments.

Color Key	
Info Complete	Canceled
Info Missing	Overlap

Calendar version 1.7.15 [2022.03.16]
 by Francis Favorini.

	Sunday 5/7/2023	Monday 5/8/2023	Tuesday 5/9/2023	Wednesday 5/10/2023	Thursday 5/11/2023	Friday 5/12/2023	Saturday 5/13/2023
8:00 AM	Add item	Add item	Add item	Add item	Add item	Add item	Add item
9:00 AM		Null.01		Null.01	Null.01		
10:00 AM		NetTMS.01 Margaret McAllister 9:30 AM-11:00 AM	Null.01 NetTMS.01 Margaret McAllister	NetTMS.01 Margaret McAllister	NetTMS.01 Margaret McAllister		
11:00 AM			Null.01 NCANDA.04 Casey Pearce 11:30 AM-1:00 PM	NOGGIN.01 Jeffrey Browndyke 11:00 AM-12:30 PM	NCANDA.04 Victoria Swaine		Null.01 Anne Baker 10:30 AM-1:00 PM 00000
12:00 PM		NoTech.01 None		Null.01	Null.01		
1:00 PM		Null.01	Null.01	NoTech.01	Null.01	NoTech.01 None	Null.01 NCANDA.04 Victoria Swaine
2:00 PM		ADRC.01 Emily Maingi	NoTech.01	Null.01			
3:00 PM		Null.01					
4:00 PM		ADRC.01 Emily Maingi					
5:00 PM		Null.01 None					
6:00 PM							
7:00 PM							
8:00 PM							
9:00 PM							



Scheduling a Test Room

To Schedule a Test Room...

- Same as Null.01 or a paid scan
- DO NOT need nulls for the test rooms
- Do not need to book within tech hours
- If no longer need slot, please cancel or delete scheduled time
- No fee for the rooms or a cancellation fee
- ****Must complete calendar entries for all test rooms used for renewal as well****

BIAC Scanner Schedule

[Experiment Info Page](#)
[MR Tech Schedule](#)

Scanner: ▼

Start Time: **Active**


Stop Time:

S	M	T
	1	2
7	8	9
14	15	16
21	22	23
28	29	30

Inactive

- BIAC1
- BIAC2
- BIAC3
- HOCK1
- MR1
- MR4
- RONC1
- VAERP
- WSPR1

Schedule just like a paid scan



BIAC Scanner Schedule

[Experiment Info Page](#)
[MR Tech Schedule](#)

Scanner:

Start Time:

Stop Time:

May

S	M	T	W	T	F	S
	1	2	3	4	5	6
7	8	9	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30	31			

Experimenters must complete all calendar entries for their experiments.

Color Key	
Info Complete	Canceled
Info Missing	Overlap

Calendar version 1.7.15 [2022.03.16]
by [Francis Favorini](#).

	Sunday 4/30/2023	Monday 5/1/2023	Tuesday 5/2/2023	Wednesday 5/3/2023	Thursday 5/4/2023	Friday 5/5/2023	Saturday 5/6/2023
	Add item	Add item	Add item	normal tech coverage until 11 11-1 - 2 techs normal tech coverage after 2	Add item	Add item	no tech
8:00 AM							
9:00 AM							
10:00 AM		NCANDA.81 Casey Pearce			NetTMS.81 Margaret McAllister 10:00 AM-11:30 AM	BrainSpine.81 Anne Baker	
11:00 AM							
12:00 PM				ALADDIN.81 Rachele Brassard 11:30 AM-1:30 PM			
1:00 PM			ALADDIN.81 Courtney Holland				
2:00 PM		UserTest.01 Kenneth Roberts					
3:00 PM							
4:00 PM					BrainSpine.81 Anne Baker		
5:00 PM		BrainSpine.81 Anne Baker		Free Slot			
6:00 PM							
7:00 PM							
8:00 PM							
9:00 PM							

*Do not need a Null – just need to schedule the time you need

Removing a scan

Cancelling and Deleting

Delete Study vs Cancel Study

- Delete
 - More than 48 hours before the scan's start time
 - No cancellation fee
 - Completely removed from calendar
- Cancel
 - Within 48 hours of the scan's start time
 - Cancellation fee (\$160 an hour)
 - Entry appears on the calendar in red

Deleting a Scan...

Study Information

biac-web.dhe.duke.edu/calendar/recdetail.asp?ScheduleID=153099

Study Details

[Experiment Info Page for Demo.01](#)

Schedule ID	153099	Study ID		Scheduled Times	5:00 PM - 5:30 PM
Date	5/2/2023 Tuesday	*Actual Start Time	5:00 PM	*Actual End Time	5:30 PM
Scanner	BIAC6	*Experiment	Demo.01	*Experimenter	None
*First Exam #		*Last Exam #		*MR Tech	
Subject		<input type="button" value="Change Subject"/>		*Subject Type	None
Consent Form Signed		Completed	<input checked="" type="checkbox"/>	*Amount Paid \$	0.00
*Problem	None			Lost Time	0 minutes
Problem Details					
Notes					

*Required fields.

...Removes a scan from calendar completely

Cancelling a scan

The screenshot shows a web browser window with the address bar displaying `biac-web.dhe.duke.edu/calendar/recdetail.asp?ScheduleID=153071`. The page title is "Study Information". The main content area is titled "Study Details" and contains a form for entering scan information. A dropdown menu is open under the "Problem" field, showing a list of options: "None", "Other_Problem", "Started_Early", "Started_Late", "Finished_Early", "Subject_Problem", "Subject_Incident", "No_Tech", "Prior_Scan_Ran_Over", "BIAC_Equip_Error", "Scanner_Error", "Scanner_Down", and "Severe_Weather". The "Other_Problem" option is currently selected. Other fields include "Schedule ID" (153071), "Date" (5/1/2023, Monday), "Scanner" (BIAC5), "*First Exam #" (01220), "Study ID" (230501_01220), "*Actual Start Time" (5:00 PM), "*Experiment" (Demo.01), "*Last Exam #" (01220), "Scheduled Times" (5:00 PM - 6:00 PM), "*Actual End Time" (6:00 PM), "*Experimenter" (None), "*MR Tech" (Jenn), "*Subject Type" (None), "*Amount Paid \$" (0.00), and "Lost Time" (0 minutes). A "Completed" checkbox is checked. At the bottom, there are buttons for "Save", "Don't Save", "Cancel Study", "Delete Study", and "Un-Cancel Study". A red message at the bottom left says "No subject".

Study Information

Study Details

[Experiment Info Page for Demo.01](#)

Schedule ID: 153071

Date: 5/1/2023 Monday

Scanner: BIAC5

*First Exam #: 01220

Subject: [Redacted]

Consent Form Signed: [Redacted]

*Problem: [Dropdown Menu]

Problem Details: [Redacted]

Notes: [Redacted]

*Required fields.

Save Don't Save Cancel Study Delete Study Un-Cancel Study

Study ID: 230501_01220

*Actual Start Time: 5:00 PM

*Experiment: Demo.01

*Last Exam #: 01220

Change Subject

Completed:

Scheduled Times: 5:00 PM - 6:00 PM

*Actual End Time: 6:00 PM

*Experimenter: None

*MR Tech: Jenn

*Subject Type: None

*Amount Paid: \$ 0.00

Lost Time: 0 minutes

No subject

Step 1) Need to select a problem from the dropdown – be as accurate as possible

Cancelling a scan

Study Information x +

biac-web.dhe.duke.edu/calendar/recdetail.asp?ScheduleID=153071

Study Details

[Experiment Info Page for Demo.01](#)

Schedule ID	153071	Study ID	230501_01220	Scheduled Times	5:00 PM - 6:00 PM
Date	5/1/2023 Monday	*Actual Start Time	5:00 PM	*Actual End Time	6:00 PM
Scanner	BIACS	*Experiment	Demo.01	*Experimenter	None
*First Exam #	01220	*Last Exam #	01220	*MR Tech	Jenn
Subject		<input type="button" value="Change Subject"/>		*Subject Type	None
Consent Form Signed		Completed	<input checked="" type="checkbox"/>	*Amount Paid \$	0.00
*Problem	Other_Problem			Lost Time	0 minutes
*Problem Details	no longer needed				
Notes					

*Required fields.

Step 2) Type in problem details

Cancelling a scan

Study Information

biac-web.dhe.duke.edu/calendar/recdetail.asp?ScheduleID=153071

Study Details

[Experiment Info Page for Demo.01](#)

Schedule ID	153071	Study ID	230501_01220	Scheduled Times	5:00 PM - 6:00 PM
Date	5/1/2023 Monday	*Actual Start Time	5:00 PM	*Actual End Time	6:00 PM
Scanner	BIAC5	*Experiment	Demo.01	*Experimenter	None
*First Exam #	01220	*Last Exam #	01220	*MR Tech	Jenn
Subject	[Redacted]	Change Subject		*Subject Type	None
Consent Form Signed	[Redacted]	Completed	<input checked="" type="checkbox"/>	*Amount Paid \$	0.00
*Problem	Other_Problem			Lost Time	0 minutes
*Problem Details	no longer needed				
Notes	<div style="border: 1px solid gray; height: 40px;"></div>				

*Required fields.

Save Don't Save **Cancel Study** Delete Study Undo Changes Un-Cancel Study

***Must click “Cancel Study” – Clicking “Save” will not cancel the scan**

****If canceled – will appear red on the scanner schedule**

Cancelling a scan – Subject Problem

- All problems involving the subject are a subject problem
 - Examples:
 - No shows or last minute cancellations
 - Implant not cleared and could not be scanned
 - Claustrophobia
 - Did not fit in the scanner
- Please be as accurate as possible and do not list subject problems as “Other Problem” or “Scanner Error”

Calendar Entries

Calendar Entries

- Who is responsible for completing calendar entries?
 - The experimenter running the study (NOT the MR Tech)
- What about if I use Mock I, Test I, DHNI, and DAVI?
 - Those entries still need to be completed!

*Incomplete calendar entries will cause your BIAC experiment to be ineligible for renewal

Calendar Items

Completed before the scan

- Actual Start and End Times
- Experimenter
- Subject Name**
- Subject Type

Completed after the scan

- Exam Number(s)
- MR Tech
- Consent Form Signed
- Problem (if applicable)
- Lost Time (if BIAC problem)

Calendar Entries: Completed Items

Study Details

[Experiment Info Page](#)

Schedule ID 121497

Study ID 190514_00000

Scheduled Times 2:00 PM - 4:00 PM

Date 5/14/2019 Tuesday

→ *Actual Start Time 2:00 PM

→ *Actual End Time 4:00 PM

Scanner BIAC5

→ *Experiment

→ *Experimenter None

→ *First Exam # 00000

*Last Exam # 00000

→ *MR Tech None

→ Subject

Change Subject

→ *Subject Type None

→ Consent Form Signed

→ *Amount Paid \$ 0.00

→ *Problem None

Completed

→ Lost Time 0 minutes

Problem Details

Notes

*Required fields.

Save Don't Save Cancel Study Delete Study Undo Changes Un-Cancel Study

Calendar Entries

- What happens if calendar entries are incomplete?
 - Incomplete calendar entries will cause your BIAC experiment to be ineligible for renewal

BIAC Scanner Schedule

[Experiment Info Page](#)
[MR Tech and Pre-Sold Slots Schedule](#)

Scanner:
 Start Time:
 Stop Time:

S	M	T	W	T	F	S
						1 2
						3 4 5 6 7 8 9
						10 11 12 13 14 15 16
						17 18 19 20 21 22 23
						24 25 26 27 28 29 30
						31

Experimenters must complete all calendar entries for their experiments.

	Sunday 5/10/2015	Monday 5/11/2015	Tuesday 5/12/2015	Wednesday 5/13/2015	Thursday 5/14/2015	Friday 5/15/2015	Saturday 5/16/2015
	Add item	No tech 8-Noon. Single tech coverage noon-10pm.	Add item	Normal tech coverage to 4pm one tech after 4pm Luke Out at 1230 pm.	One tech 2-6, no tech after 6.	Add item	Add item
8:00 AM			UserTest.01	QA.04			
9:00 AM			ConUp.01 Emma Dowd 8:30 AM-10:00 AM 19765			SAT.01 Hanna Oh 8:30 AM-10:00 AM 19775	
10:00 AM			UserTest.01 Jessica Stanek	UserTest.01 Jessica Stanek	SAT.01 Hanna Oh 9:30 AM-11:00 AM 19773	UserTest.01	QA.04
11:00 AM				ConUp.01 Emma Dowd 10:30 AM-12:00 PM 19771		Predator.03 Chelsea Swanson 11:00 AM-12:30 PM Canceled	
12:00 PM		QA.04				QA.04	
1:00 PM		DNS.01 Devin Ulrich					
2:00 PM		Predator.03 Chelsea Swanson 2:00 PM-3:30 PM Canceled		DECIDE.01 Puja Patel 2:30 PM-4:00 PM 19772		VisAtten.02 Catherine Tallman 2:30 PM-4:00 PM 19778	
3:00 PM		Null.01 Kayla Hernandez 3:30 PM-5:30 PM 00000	UserTest.01	DECIDE.01 Puja Patel		Null.01 Kayla Hernandez 4:30 PM-6:00 PM 00000	
4:00 PM			UserTest.01	Null.01 Jeffrey Browndyke 5:00 PM-6:30 PM 00000	Null.01 Hannah Hayward 4:00 PM-6:00 PM 00000		
5:00 PM		ConUp.01 Emma Dowd 5:30 PM-7:00 PM 19764	PAC.01 Cory Campbell				
6:00 PM			QA.04				
7:00 PM							
8:00 PM							
9:00 PM							

Color Key
Info Complete Canceled
Info Missing Overlap

Calendar version 1.7.8 [2018.09.17]
 by [Francis Favorini](#).

Calendar Entries: Problems

- Most Common:
 - Subject Problem
- Are not charged for:
 - No Tech
 - Scanner Error
 - BIAC Equipment Error
 - Scanner Down
 - Severe Weather

Schedule ID: 153071
Date: 5/1/2023 Monday
Scanner: BIAC5
*First Exam #: 01220
Subject: [Redacted]
Consent Form Signed: [Redacted]
*Problem: None
Problem Details: [Dropdown Menu]
Notes: [Redacted]

Buttons: Cancel Study, Delete Study, Un-Cancel Study

Claiming Lost Time

- Lost Time should only be claimed if you are unable to use all or any part of your data due to a BIAC-related problem
- BIAC will enter lost time for the following situations:
 - No Tech Available*
 - BIAC Scanner Down
 - Severe Weather (Duke Policy Activated)

Scheduled Times	5:00 PM - 6:00 PM
*Actual End Time	6:00 PM ▼
*Experimenter	None ▼
*MR Tech	Jenn ▼
*Subject Type	Phantom ▼
*Amount Paid \$	0.00
Lost Time	0 minutes

*****If your scan encountered a BIAC problem, please be sure to submit a trouble ticket ASAP so that the problem can be investigated and avoided for future scans*****

Trouble Tickets

- Trouble tickets can be submitted by clicking on the “Services” tab on the BIAC Homepage or by visiting www.biac.duke.edu/services/

The screenshot shows the BIAC website interface. At the top, there is a navigation bar with links for "Calendar", "Directory", and "Volunteer". Below this is a dark blue header with the "biac" logo and the text "Duke-UNC Brain Imaging and Analysis Center". A search bar is located on the right side of the header. Below the header is a navigation menu with options: "About Us", "Research", "Education", "People", "Facilities", and "Services". The "Services" option is selected, and a breadcrumb trail shows "Home > Services". The main content area features four service tiles: "Scanner Calendar" (with a calendar icon), "Trouble Ticket" (with a yellow ticket icon labeled "TROUBLE"), "Exam Tracker" (with a "LOG" icon showing "01/01 MESSAGE"), and "Status Monitor" (with a database icon). The "Trouble Ticket" tile is circled in red.

Common Scheduling/Calendar Errors

- Not knowing who is responsible for what
- Scheduled without subject
- Experiment has no hours left
- Experiment/IRB protocol has expired
- Improper Null.01 use
- Improper UserTest.01 use
 - This is NOT a placeholder for an experiment
- Incomplete calendar entries
 - Missing exam number
 - MR Tech not listed
 - Lost time not claimed properly
 - MockI/TestI/DHNI/DAVI entries are not completed
 - Consent form

Policy on Excessive Scheduling Errors

- If a user/scheduler fails to schedule Null.01 during single tech coverage or schedules all/part of their session outside of the normal or posted tech hours the following will happen:
 - **1st Offense:** Scheduler gets notified about the policy and New User Orientation.
 - **2nd Offense:** MR Operations Director, Dr. Todd Harshbarger, notifies the Scheduler and the Responsible PI about the severity of the situation, the second offense, and that on the third offense scheduling rights will be taken away.
 - **3rd Offense:** BIAC Director, Dr. Allen Song, notifies the Scheduler and the Responsible PI that the third offense has occurred and that the scheduler's scheduling rights have been removed.

Protocol Planning

UserTest.01

- UserTest.01 = free testing time
- UserTest.01 = no MRI data, but you can...
 - Set up your protocol
 - Test scanner triggering
 - Time your task
 - See what your stimuli looks like inside the scanner
- UserTest.01 can be scheduled <24 hours in advance
 - Paid study CAN schedule over a User Test!
- We highly recommend scheduling a UserTest.01 prior to your first participant!

Scan Session: Planning your Time

- A scheduled session must include:
 - (1) **Protocol time**: the time it takes to run a protocol
 - (2) **Breaks**: if there is a short resting period between runs
 - (3) **Set-up time**: the time it takes to review the screening forms and getting the participant set up
 - (4) **Exit time**: the time it takes to get a participant out of the scanner
- Typically allot 15 min beyond your protocol time

Data Quality and Movement

- When screening participants before the MRI, incorporating a mock scan can help identify if participants will be able to remain still for the duration of the scan
- Be observant of medical conditions that may cause greater than normal movement (ex. Bronchitis in the winter months or a common cold)
- Experimenters at the scanner should be observing data acquisition and noting if there are any artifacts
 - If participant movement – remind subject to be still, or if it continues, end scan early and get a partial cancellation fee back
 - If artifact from MRI – immediately submit a trouble ticket for review from BIAC Help

Data Quality and Movement

- Utilize free UserTest sessions on calendar to run through protocol multiple times before beginning data collection
- Consult with BIAC to ensure scanning protocol is correct before booking scanner time
- Data should be reviewed for QA **as soon as possible** after collection
- Utilize the online trouble ticket system to alert BIAC to potential issues with scanning or reach out to Dr. Todd Harshbarger
 - BIAC is not responsible for artifacts related to participant motion

Participant Screening and MR Safety

Pre-Scan Screening

- All subjects must be either phone screened or interviewed for MRI safety in person before being scheduled for a scan
 - Highly recommended to administer official BIAC MRI Safety Screening Form

- Any questions or concerns about MRI safety found in the initial screening should be addressed with Todd and the MRI Techs at least 48 hours before the scan

Pre-Scan Screening

MRI Techs typically are meeting participants for the first time the day of the scan and are a last line of defense

- Responsibility to determine initial MRI safety and participant quality is up to individual lab faculty and staff (unless using BIAC volunteer registry)
- Thorough review of medical history
 1. May require obtaining medical record authorization to review make and model of any implanted metal or part of medical history
 2. Emphasis on reviewing and researching any past surgeries as many participants may not know they have implanted metal after certain procedures (i.e. wire sutures or a screw/pin)
 3. Even if implants are deemed safe for scanning, they may cause artifacts

Examples of Pre-Screening Resources

■ Phone Screening Script

BIAC Phone Screening

Part I: Initial Contact

Good morning/afternoon, my name is Marissa Bruce and I'm calling from the Duke Brain Imaging and Analysis center. Are you free now to hear more about our studies?

If yes,

Great. Would you be free to do the phone screening now?

If no,

Is there a better time for you to do the phone screening?

Part II: BIAC Research Summary

If yes,

Our research uses functional magnetic resonance imaging to take pictures of your brain while you are doing a psychology experiment. We study the areas of the brain associated with things like memory, attention, and perception. We have a number of studies going on at any one time, and our studies are conducted at the MRI center at Duke Hospital.

If you participate in our studies, you will come to the MRI center and lie down on your back in a MRI machine (which is like a small open-ended tunnel).

The machine uses radio waves to take the pictures of your brain. The radio waves and magnetic fields used in MRI are not themselves harmful, but they do pose a risk in some cases, such as if you have a pacemaker.

Most studies take between 1-2 hours, and pay about \$20 per hour, but the duration and payment depends on the particular experiment. Are you still interested in participating?

If yes,

In order to participate and be entered into our database (so that you can be contacted for our studies), the first step is a 5-10 minute phone screening and the second step is to come into the Duke Children's Hospital for a one hour compensated in-person screening. You will receive \$20 and a parking pass for the hospital garage. During the screening, we will collect a health history, have you fill out some questionnaires to see if you are eligible for our studies, and conduct an MRI simulation. Would you be interested in this screening?

If no,

Thank you for taking the time to reach out to us and have a great day!

Part II: Screening Questions

If yes,

Before you can participate in the in-person screening, I need to go over a set of questions with you. If you do not want to answer a question, you do not have to answer it, but then you will not be able to participate.

1. Have you done any previous MRI studies or ever been in a MRI scanner?

"Yes". Ask the subject when and where the session was conducted. Note: prior participation in MRI does not ensure safety!

"No." Be especially careful with any subject who has not been in the scanner.

Revised: 03/04/2019

2. Do you have any of the following metal objects in your body?

a. Implants: pacemakers, cochlear/ear, insulin pump, etc.

b. Instruments: metal working, bullet/shrapnel, etc.

"Yes". If they have any implanted metal, then they should not participate in the study or go into the scanner room without approval from Dr. Song. Remember, the scanner room itself can be a dangerous environment for subjects with implanted metal (e.g., pacemakers), so it is important to screen these subjects ahead of time. *[STOP SCREENING – CANNOT PARTICIPATE]*

3. Do you have any of the following metal objects in your body?

a. Surgical: pins in bones, clamps, valves, catheters, etc.

"Yes". Inquire about the type of implant and the type of metal. If they have any implanted metal that is not MR-compatible then they should not participate in the study or go into the scanner room without approval from Dr. Song. Remember, the scanner room itself can be a dangerous environment for subjects with implanted metal (e.g., pacemakers), so it is important to screen these subjects ahead of time.

4. Have you ever worked with metal (e.g. grinding, fabricating)

"Yes". Ask the subject if he or she was ever around metal shavings. If "Yes". *[STOP SCREENING – CANNOT PARTICIPATE]*

5. Have you ever had an injury to the eye involving metal or invasive eye surgery?

"Yes". *[STOP SCREENING – CANNOT PARTICIPATE]*

6. Are you pregnant or think you may be pregnant? (women only)

"Yes". *[STOP SCREENING – CANNOT PARTICIPATE]*

7. Do you have any permanent makeup or tattoos?

"Yes". *[STOP SCREENING – CANNOT PARTICIPATE, tattoos below the neck are ok]*

8. Do you have any body piercing that cannot be removed?

"Yes". *[STOP SCREENING – CANNOT PARTICIPATE]*

"No. I have piercings but they are removable." Tell the subject that they should leave all piercings and jewelry at home, or in lockers outside the MR center.

The following questions may preclude the subject's participation in some studies, or may suggest some potential problems:

9. Do you have dental braces, dental implants, or a permanent retainer?

"Yes. I have permanent metal dental work." We currently do not run subjects with permanent dental work because of the potential for an imaging artifact. In addition, some permanent retainers are held in place by magnets.

"Yes, but it is removable." All removable metal dental work should be removed before entering the MR room.

"Are fillings OK?" Gold and silver fillings are OK for entering the scanner room.

10. Do you feel anxious in small spaces? Would you classify yourself as being claustrophobic? (*emphasize*)

"Yes". Explain the MR environment to them, and describe the study. Tell them that they will in the scanner bore for at least 90 minutes, but will be able to talk with the experimenters in the next room. Let the subject make a decision about whether to participate.

11. Do you use any prescription medications, for example, anti-depressants, anti-anxiety or any hypertension medications?

Because of the potential effects of some medications upon brain neuronal activity and blood flow, we also ask subjects about the drugs they might be taking

Revised: 03/04/2019

Part III: Result of Phone Screening

If they qualify,

That is all of the questions. At this point, if you are still interested, we can schedule you for a visit to our MRI facility at Duke for the in-person screening.

If they do NOT qualify,

Unfortunately, you are not eligible to participate in our MRI studies. Because these studies are for research and of no medical benefit to you we must take additional precautions to ensure your safety. If you have any further questions or concerns, you can contact our director, Dr. Allen Song at 681-9337.

Part IV: Voicemail

Hello,

This is Marissa Bruce from the Brain Imaging and Analysis Center at Duke. I am calling to conduct a quick 5 minute phone screening to see if you are eligible for our studies. If you would like, please give us a call back at 919.681.9344 any time Monday through Friday between 8am-4:30pm. If you prefer, please send an email to volunteer at b-i-a-c dot duke dot edu. That is b-o-i-l-a-n-t-e-e-r at b-i-a-c dot duke dot edu. Thank you for your time!

Revised: 03/04/2019

Examples of Pre-Screening Resources

■ BIAC MRI Safety Screening Form

Protocol: _____ Exam Number: _____ Date: _____

**Duke-UNC Brain Imaging and Analysis Center:
MRI Safety Screening**

All individuals entering the MRI suite must fill out this information to the best of their knowledge. Any potential contraindications must be reviewed with the individual's medical record and the BIAC MR Safety Committee before being cleared to enter the scanner bore.

Part I: For all individuals entering the scanner room

Name _____ Birthdate _____
Last name First name M.I.

Address _____ City _____

State _____ Zip Code _____ Phone (H)(____) _____ (W)(____) _____ (C)(____) _____

1. Have you ever had an injury to the eye involving a metallic object (e.g. metallic slivers, shavings, foreign body)? No Yes
 If yes, please describe: _____

2. Have you ever worked with metal (grinding, fabricating, etc.)? No Yes
 If yes, please describe: _____

3. Have you ever had surgery (including eye surgery)? No Yes
 If yes, please describe: _____

4. Have you had any previous MRI studies or been in a MR scanner? No Yes
 If yes, please list (most recent first): Body part _____ Date _____ Facility _____
 If yes, did you have any problems? _____

Before you may enter the scanner room, you must remove all metallic objects.

<input type="checkbox"/> All contents of pockets, including back pockets	<input type="checkbox"/> Shoes that contain any metal (e.g., steel tipped)
<input type="checkbox"/> Wrist watch, any bracelets	<input type="checkbox"/> Hearing aids or other electronic devices
<input type="checkbox"/> Hair pins, clips, weaves, fasteners	<input type="checkbox"/> Pagers, cell phones, PDAs
<input type="checkbox"/> Pins or badges on shirt	<input type="checkbox"/> Dentures or removable retainer
<input type="checkbox"/> Belt with metal (e.g., buckle)	<input type="checkbox"/> Necklaces, chains

Part II: For all individuals entering the scanner bore

1. Are you claustrophobic? No Yes

2. Do you have an IUD or diaphragm containing metal? No Yes

3. Are you pregnant, experiencing late menstrual period, or undergoing fertility treatment? No Yes

4. Do you currently have a fever or other acute illness? No Yes

5. Please list any surgeries or other invasive medical procedures in as much detail as possible:

Revised: 7/28/2017 Page 1 of 3

Protocol: _____ Exam Number: _____ Date: _____

6. Are you currently taking or have you recently taken any medication? No Yes
 If yes, please list _____

7. Do you have anemia or any diseases that affect your blood? No Yes
 If yes, please describe _____

8. Do you have a history of stroke, seizures, brain tumor, head trauma, or other neurological disorder? No Yes
 If yes, please describe _____

9. Do you wear glasses or contact lenses? No Yes
 If yes, please specify prescription (if known) _____

10. Do you have a breathing disorder (e.g., asthma, apnea), heart condition, or movement disorder? No Yes

Height _____ Weight _____ Handedness _____

WARNING: Certain implants, devices, or objects may be hazardous to you and/or may interfere with the MR procedure (i.e., MRI, MR angiography, functional MRI, MR spectroscopy). Do not enter the MR system room or MR environment if you have any questions or concerns regarding an implant, device, or on object. Consult the MRI Technologist or Radiologist BEFORE entering the MR system room. The MR magnet is ALWAYS on.

Please mark on the figure(s) below the location of any implant or metal inside of or on your body.

Revised: 7/28/2017 Page 2 of 3

Protocol: _____ Exam Number: _____ Date: _____

Please indicate if you have any of the following:

<input type="checkbox"/> Yes <input type="checkbox"/> No	Aneurysm clip(s)
<input type="checkbox"/> Yes <input type="checkbox"/> No	Cardiac pacemaker
<input type="checkbox"/> Yes <input type="checkbox"/> No	Implanted cardioverter defibrillator (ICD)
<input type="checkbox"/> Yes <input type="checkbox"/> No	Electronic implant or device
<input type="checkbox"/> Yes <input type="checkbox"/> No	Magnetically-activated implant or device
<input type="checkbox"/> Yes <input type="checkbox"/> No	Neurostimulation system
<input type="checkbox"/> Yes <input type="checkbox"/> No	Spinal cord stimulator
<input type="checkbox"/> Yes <input type="checkbox"/> No	Internal electrodes or wires
<input type="checkbox"/> Yes <input type="checkbox"/> No	Bone growth/bone fusion stimulator
<input type="checkbox"/> Yes <input type="checkbox"/> No	Cochlear, otologic, or other ear implant
<input type="checkbox"/> Yes <input type="checkbox"/> No	Insulin or infusion pump
<input type="checkbox"/> Yes <input type="checkbox"/> No	Implanted drug infusion device
<input type="checkbox"/> Yes <input type="checkbox"/> No	Any type of prosthesis (eye, penile, etc.)
<input type="checkbox"/> Yes <input type="checkbox"/> No	Heart valve prosthesis
<input type="checkbox"/> Yes <input type="checkbox"/> No	Eyelid spring or wire
<input type="checkbox"/> Yes <input type="checkbox"/> No	Artificial or prosthetic limb
<input type="checkbox"/> Yes <input type="checkbox"/> No	Metallic stent, filter, or coil
<input type="checkbox"/> Yes <input type="checkbox"/> No	Shunt (spinal or intraventricular)
<input type="checkbox"/> Yes <input type="checkbox"/> No	Vascular access port and/or catheter
<input type="checkbox"/> Yes <input type="checkbox"/> No	Radiation seeds or implants
<input type="checkbox"/> Yes <input type="checkbox"/> No	Medication patch (Nicotine, Nitroglycerine)
<input type="checkbox"/> Yes <input type="checkbox"/> No	Any metallic fragment or foreign body
<input type="checkbox"/> Yes <input type="checkbox"/> No	Wire mesh implant
<input type="checkbox"/> Yes <input type="checkbox"/> No	Tissue expander (i.e. breast)
<input type="checkbox"/> Yes <input type="checkbox"/> No	Surgical staples, clips, or metallic sutures
<input type="checkbox"/> Yes <input type="checkbox"/> No	Joint replacement (hip, knee, etc.)
<input type="checkbox"/> Yes <input type="checkbox"/> No	Bone/joint pin, screw, nail, wire, plate, etc.
<input type="checkbox"/> Yes <input type="checkbox"/> No	Dentures or partial plates
<input type="checkbox"/> Yes <input type="checkbox"/> No	Tattoo or permanent makeup
<input type="checkbox"/> Yes <input type="checkbox"/> No	Body piercing or jewelry
<input type="checkbox"/> Yes <input type="checkbox"/> No	Hearing aid (remove before entering MRI)
<input type="checkbox"/> Yes <input type="checkbox"/> No	Other implant

IMPORTANT INSTRUCTIONS

Before entering the MR environment or MR system room, you must remove all metallic objects including hearing aids, dentures, partial plates, keys, cell phone, eyeglasses, beeper, hair pins, barrettes, jewelry, body piercing jewelry, watch, safety pins, paperclips, money clip, credit cards, bank cards, magnetic strip cards, coins, pens, pocket knife, nail clipper, tools, clothing, with metal fasteners, and clothing with metallic threads. You will be asked to wear ear plugs to protect your hearing during the scan.

Please consult the MRI Technologist or Radiologist if you have any question or concern BEFORE you enter the MR system room.

I attest that the above information is correct to the best of my knowledge. I read and understand the contents of this form and had the opportunity to ask questions regarding the information on this form and regarding the MR procedure that I am about to undergo.


Signature of Person Completing Form _____ Signature of Person Screening Subject/Patient _____ Date _____

Form Completed by: Self Parent/guardian Other relative Physician

Revised: 7/28/2017 Page 3 of 3

MR Safety

MR Safety: Zones



- **Zone I**
 - All areas accessible to the general public
- **Zone II**
 - Areas between the public-accessible Zone I and restricted areas of Zone III and IV
- **Zone III**
 - Area where free access by unscreened non-MR personnel or ferromagnetic objects or equipment can result in serious injury or death
- **Zone 4**
 - MR scanner room


MR Safety

- 24 hrs/day
- 365 days/year
- The magnetic field is on even when the magnet is **NOT** in use.
- Even if the electrical portion of the machine is off, the magnetic field is still up
- The MRI technologist has final say in who may be allowed to go into the operation and/or scanner room.



MR Safety

- The force of the field is measured in Tesla (T)
 - Our scanners are measured at 3.0 Tesla (T)
- The force of the field is greatest at the center of the magnet.
 - The magnetic FORCE INCREASES as you move closer to the center of the magnet
- What is not attracted to the magnet outside the fringe field can be attracted to the magnet as you move closer to the center of the machine



Participant Safety is our #1 Concern!

MR Safety: Medical Implants

- Medical Implants
 - Review Dr. Frank Shellock's MRI contraindication database: www.mrisafety.com
 - At least 3 days before your exam, email the device model and manufacturer of any medical implants to all 3 techs at BIAC-Schedule@duke.edu
 - They will reply via email regarding eligibility
 - It is the Experimenter's responsibility to get approval prior to the scan!

Magnetic Resonance Imaging
Tips for Scanning Patients with Implants

• Follow your site's process for screening the patient
• Identify the manufacturer and model of any implanted devices
• Locate the MRI safety information in the device manufacturer's labeling

Look for one of these icons:

BEFORE

- MR** MR Safe. Patients with MR Safe devices have no scanning restrictions.
- MR** MR Conditional. For patients with MR Conditional devices, implant conditions should be matched with the MR system information.
 - Consult your MR system manual for MR system information
 - Ensure that the MR system meets all conditions provided in the MR Conditional labeling
 - If conditions are not met, the patient should not be scanned
- MR** MR Unsafe. Patients with MR Unsafe devices should not be scanned. Assume any unidentified implant is MR Unsafe.

• Document device information in the medical record
• Consult a physician for any risk/benefit decisions
For MR Conditional devices:
• Follow all pre-scan conditions, such as special programming modes

DURING

- For MR Conditional devices, follow all scan conditions such as specific absorption rate (SAR) restrictions or patient positioning instructions
- Monitor the patient at all times

AFTER

- Assess the patient for discomfort or injuries
- Follow any post-scan conditions, such as device checks or programming


FDA U.S. FOOD & DRUG ADMINISTRATION
SMRT A WORLD OF KNOWLEDGE PROFESSIONALS

MR Safety: Pregnancy Tests

- Although BIAC no longer requires pregnancy tests, many studies still have this requirement in their IRB procedures.
- Female participants (12-55 years) must receive a urine immediately prior to each MRI Scan
- Very limited exceptions:
 - Females who have not yet menstruated
 - Post-menopausal females
 - Females who have had a hysterectomy
- Supplies are available at the scanners
- Experimenter should facilitate urine collection
- MR Tech will perform pregnancy test



Pregnancy Test Stations



Pregnancy Tests Update

- Since February, 2023, BIAC no longer requires pregnancy tests, so it is now up to each study to keep or get rid of this requirement
- If the study no longer wants to require pregnancy tests, then they must submit an amendment to their IRB
- If the study would like to continue requiring pregnancy tests, BIAC will still provide kits and testing

What to expect the day of the
scan...

Day of Scan Procedure

Step 1: Escort participant to MRI hallway and ask them to place all belongings into lockers

Step 2: Ask participant to remove all metal objects from their person and place in locker (jewelry, bobby pins, etc.)

Step 3: Ask participant to pat and empty any and all pockets

All participants should be MRI ready before entering the console room

Step 4: If scanning a woman of child bearing potential, conduct the urine pregnancy test

Step 5: Walk participant to console room and have them fill out the *BIAC MRI Safety Sheet*

Step 6: Ask participant to re-check for any metal on their person or in their pockets

Step 7 : MRI Tech will review form and begin their scanning procedure

Scan Session:Visitors

- Visitors and family members are never allowed in the console room
- What are your options?
 - Schedule time in Test I, DHNI, or Mock I
 - Parent/guardian can wait in the Children's Hospital Lobby or the Main Entrance Lobby
- **CANNOT** leave family/visitors:
 - In the hallway
 - Console room



Scan Session: Hospital Courtesies

- At the hospital, patients and their transporters have the right of way
- Please be sure to wait if patients are in the middle of the MR hallway
- ID Badge must be visibly worn at all times

Scan Session: Experimenter Courtesies

- Please be sure to wait until the scan before yours is finished, and the participants have exited the console room before entering
- What are your options?
 - Schedule TestI/DHNI/MockI* and wait in those rooms
- Please do not wait in the hallways
- Scans are not allowed to run over the scheduled time

Scan Session: Experimenter Courtesies

- Do NOT train/consent/test participants in the hallways
 - This includes the table by BIAC 5
- If it is necessary to talk with a member from another lab, pull them aside and do not talk to them in front of their participants
- Additionally, a participant may be from a sensitive or vulnerable population and require the undivided attention of the lab team.

BIAC Resources

BIAC Resources: Stimulus Controls

- We have many peripherals available for you to use!

- Details instructions can be found here:
<https://wiki.biac.duke.edu/biac:experimentalcontrol>

BIAC Stimulus Control System Manual

Stimulus Control System Layout

- Organization diagram
- Hardware connected to each PC

Video display

- Using the video projector
- Using the DVD player
- Using the Video switch
- Using the Monitors

Audio

- Using the audio headphones

EyeTracker

- Eyetracking in the MRI scanner

Galvanic skin resistance (GSR)

- Recording GSR in the MRI scanner

Button Boxes

- Using the 8-button box
- Using the joystick button box

Head Tracking Device

- Using the head tracker

Troubleshooting

- Troubleshooting


Table of Contents

- BIAC Stimulus Control System Manual
- Stimulus Control System Layout
- Video display
- Audio
- EyeTracker
- Galvanic skin resistance (GSR)
- Button Boxes
- Head Tracking Device
- Troubleshooting

BIAC Resources: Participant Recruitment

- BIAC maintains an active participant registry (~100 members)
- Healthy, “Normal” participants
- Criteria:
 - 18-80 years of age
 - Right-handed
 - Earned at least a high school diploma
 - English speaker (native/fluent)
 - No history of brain injury or disease
 - No anti-depressants, anti-anxiety, or hypertension medication use

****Currently not recruiting, but can still use**




Participants Needed
Duke Brain Imaging and Analysis Center

Healthy participants are needed for research studies of visual and hearing function using magnetic resonance imaging (MRI) at Duke University Medical Center. Most studies last between **1-2 hours** and **participants will be paid** for their time.

PARTICIPANT REQUIREMENTS

- 18-80 years of age
- Right-Handed
- Earned at least a high school diploma
- Fluent in English
- No history of brain injury or disease

 **DukeMedicine** Learn more at biac.duke.edu
919-681-9344 | volunteer@biac.duke.edu

BIAC Resources: Using BIAC Participant Registry

If needed, submit IRB amendment to add BIAC Participant Registry Protocol (Pro00010672) to your protocol

- Note: protocol must be approved by the BIAC Scientific Review Committee prior to beginning recruitment

If needed, provide study title, PI name(s), and IRB number to the BIAC Research Coordinator for IRB submission

Once both protocols have been approved, submit either a “Behavioral Study Request” form or a “MRI Study Request” form to the BIAC Research Coordinator

Please confirm that you have requested Mock I/Test I/ DHNI/DAVI on your BIAC Research Proposal if you need a room for consenting, debriefing, or behavioral tasks

BIAC’s Research Coordinator will use the provided information to contact eligible participants from our Participant Registry

BIAC Resources: Using BIAC Participant Registry

- Please be aware that we schedule and recruit for many BIAC studies. Please provide enough time for IRB amendment approval and/or contacting participants if you need subjects.
- Please try to be flexible with your schedule, our scanners have extended hours Monday - Wednesday nights, plus availability on Thursday nights/Saturdays alternate.
- Keep in mind that many of our participants are students, be mindful of the academic calendar (fall break, exams weeks, major holidays, spring break, etc.)

BIAC Resources: Participant Feedback

- Please let us know if a participant does not show for a scan as they will no longer be contacted for our studies
- Please let us know if a participant performs extremely poorly on a task, is unusually unmotivated, or falls asleep during the scan
 - We have a two-strike policy with participants for instances like these, so the more feedback—the better our recruitment can be
- Please let us know if the participant has trouble laying still
 - Dental work may produce image artifacts

BIAC Billing System

BIAC Billing System

- Each month, we generate invoices for each experiment through the BIAC Scanner Schedule System
 - Invoices will be emailed on the 1st business day of the following month
- Calendar entries must be fully documented by the end of the last day of every month. Once the invoices are generated, changes cannot be made.
 - This includes making any lost time claims and/or cancellations

BIAC Billing System

- Invoices are sent to the following individuals that were listed on the initial BIAC proposal form received by the Scientific Review Committee
- These individuals should review the invoices upon receipt and contact billing@biac.duke.edu to confirm the charges
- Payments will be processed within 2-4 days after the invoices are sent out
 - We will charge the last cost object (fund code) that we have on file for the study. If the funding source has changed, please notify BIAC as soon as you have the new fund code.
- Updating calendar entries regularly will avoid the need for corrections on the invoices.

BIAC Billing System

BIAC Scanner Time Invoice

Please remit payment to:

Center Attn: Billing Campus Box 3918 Durham, NC 27710	Brain Imaging and Analysis 919-681-7033 Fax billing@biac.duke.edu	919-681-9337
--	---	--------------

Department **BIAC**

Responsible Investigator **Allen Song** allen.song@duke.edu
Billing Contact **Elena Musty** elena.musty@duke.edu

Experiment [REDACTED] **Cost Object Investigato** [REDACTED]

Scanner	Date	Start	End	Exam(s)	Canceled	Lost	TimeOverlap	Amount	Paid to	Amount	
								Credit	Billed	Date	Due

Current Charges

BIAC6	15-Apr-2019	2:30 PM	3:00 PM	00606	No	0 min.	0 min.	\$280.00	\$0.00	\$280.00
BIAC6	25-Apr-2019	3:30 PM	4:30 PM	00645	No	0 min.	0 min.	\$560.00	\$0.00	\$560.00

Current charges through 4/30/2019 (2 studies)	\$840.00
Past Due charges (0 studies)	\$0.00
Development funds contributed by BIAC	\$0.00
Total for [REDACTED] (2 studies)	\$840.00

Current charges through 4/30/2019 (2 studies)	\$840.00
Past Due charges (0 studies)	\$0.00
Development funds contributed by BIAC	\$0.00
Total for Allen Song (2 studies)	\$840.00

BIAC Billing System

- Rates Effective July 1st, 2023:
 - Rate for Scans = \$640 per hour
 - Cancellation Fee Rate = \$160 per hour
 - Scans cancelled with less than 48 hours notice will incur the cancellation fee

BIAC Billing System: Helpful Definitions

- **Lost Time:** This is a credit for time/data lost due to a BIAC-related problem. The experimenter should claim this on the BIAC Scanner Schedule
- **Overlap Credit:** This is a credit granted when a paying study takes over a canceled slot and completes the scan
 - Non-Paying Studies:
 - UserTest
 - SysDevelop
 - SysService
 - QA
 - Demo

BIAC Billing System: Cancellations

- Canceled scans incur a fee at the rate of \$160/hr. (effective 7/1/2023).
 - Please note that this charge will only apply to those sessions which are canceled <48 hours of the session start time
- Sessions canceled with >48 hours notice, or due to scanner down time, will have no associated charges. If a BIAC MR Tech becomes unavailable for your scan, we will not charge you for the time.
- The cancellation fee can only be applied to the entire session
 - You will be charged the full rate for your scheduled time, even if you only use a portion of it.
- The earlier you cancel your study on the calendar, the more time another paying study has to take over the slot

FAQs

FAQs

- Q: I am at the scanner and my session ended and I need more time. What do I do?

FAQs

- Q: I am at the scanner and my session ended and I need more time. What do I do?
 - A: If the MR Tech and scanner are available (no other session booked), you may continue your scan by adding additional time to the calendar. If you do not add it, the MR Tech will and you will be billed appropriately. You can only add time in 30 minute increments. (i.e. If you schedule 30 minutes but need 45, you will pay for 1 hour. The time you schedule on the scanner should be adequate to cover time you enter the console room, until the time you leave. This includes the time for the participant to enter and exit the scanner.

FAQs

- Q: I no longer need my full time slot, can I cancel a portion of my scheduled scan time for a reduced fee?

FAQs

- Q: I no longer need my full time slot, can I cancel a portion of my scheduled scan time for a reduced fee?
 - A: You can cancel the portion on the calendar, however you will not receive a reduced fee unless a paying study takes over your time slot (overlap credit will apply). If you scheduled time and run your study, you will be charged for the entire scheduled slot regardless of how much time you use.

FAQs


- Q: My participant has an implant/device/object but says they have had a MRI scan since then. Do I still need to get the MR Techs to clear it before their scan? Even if the scan was done at Duke?

FAQs

- Q: My participant has an implant/device/object but says they have had a MRI scan since then. Do I still need to get the MR Techs to clear it before their scan? Even if the scan was done at Duke?
 - A: YES! You will still need to have the implant/device/object cleared by the MR Techs at least 3 days prior to the scan.
 - A: YES! Both 1.5T and 3T scanners are used clinically, so implants/devices/objects that were at one time cleared for 1.5T have not all been cleared for 3T.

Next Steps...

Next Steps...

1. **Get a DHE account**
2. Send NetID and email to biac-newuser@duke.edu
3. Complete the BIAC User Profile survey
4. Complete the Level I Access survey
 - Date of scheduled New User Orientation
 - Date completed MRI Safety Quiz
 - Uploaded MRI Level I Access Signature Sheet
 - Upload MRI Subject Screening form
5. **Complete the New User Orientation session** 
6. Join the BIAC User Listserv

Forms and Safety Quiz

- Found at <https://www.biac.duke.edu/research/safety>
- Uploaded to the Level I Access survey

Safety

If you need access to the MRI suite, you must review the MRI safety materials and pass the safety quiz to demonstrate your understanding of the material covered. You must also complete, sign and date the MRI Screening Form and Level 1 MRI Access Signature Sheet and return them to the [BIAC Main Office](#).

- [SOP101: MRI Access](#)
- [SOP111: BIAC Safety Training](#)
- [MRI Safety Tutorial](#)
- [MRI Safety Video](#)
- [MRI Safety Quiz](#)
- [MRI Subject Screening Form \(PDF\)](#)
- [MRI Level 1 Access Signature Sheet \(PDF\)](#)

BIAC User Listserv

- Join to stay up to date and receive any notices, policy changes, etc.
- To join:
 - Go to <https://lists.duke.edu/sympa>
 - Log in with Net ID and click “Search for lists”
 - Search “biac” and click the group biac-users@duke.edu
 - On the left side of the page, click “Subscribe”

biac-users - For anyone who uses the services of the Brain Imaging & Analysis Ctr.

biac-users@duke.edu

Subject: For anyone who uses the services of the Brain Imaging & Analysis Ctr.

Description: This list is for researchers who use the BIAC facilities. Operational announcements, service interruption notices, talk announcements, and other important messages that affect the general BIAC community will be sent through this listserv.

Questions?

Thank you for attending the
BIAC New User Orientation!

